


## WHERE DO THEY COME FROM ?

HYPERTENSION,  
STENOCARDIA,  
ASTHMA  
and many others..

Low level of  $CO_2$  in blood leads to constriction of arteriols and bronchioles..

THE ROOT OF CHRONIC DISEASES  
IS LACK OF  $CO_2$  IN ARTERIAL BLOOD

## New Millennium Anti-Aging Decision

# COSMIC HEALTH

The main part of cosmic health – capnicator – serves to form the active breathing air composition, which has less oxygen but more carbon dioxide (this is regulated) than normal air. Breathing through the capnicator for 30 minutes per day for the time long about 6 months, without any tension, brings to decreasing of breathing level and, consequently, to increasing of  $CO_2$  in the blood to the optimum level 6.5%.

The second device – capnometr – is intended for measurement of the breath level. Using this device periodically (1-2 per month) you may observe the process of normalizing of your breathing.

- The usage of the cosmic health lets by 2-3 times decrease the level of breathing. The level of  $CO_2$  will increase to norm 6.5%. The tonus of microvessels will be such as at 25 years old. It is real rejuvenation!
- In the course of training it is gradually lowered the peripheral resistance of vessels and, consequently, the superfluous blood pressure.
- The normalization of blood current restore normal supplying of cells and tissues, and consequent chronic diseases.

WWW.SALUS.LT

E-mail [info@salus.lt](mailto:info@salus.lt)

More information you can get on site  
Please, send your questions on e-mail  
Address of the company : 62a, En

Tel./faks. 8-5 270305808

Mob. tel. +370 60689669

3051

## How it works

"Good breathing, like good nutrition, must meet the organism's needs".

K.P. Buteyko

Russian doctor Butejko was the first, paid attention on role  $\text{CO}_2$  with what quintessence of opening Butejko? "Breathe less", is its advice. It has drawn attention to importance of carbonic gas ( $\text{CO}_2$ ) for a metabolism of an organism and as shortage in structure of arterial blood can cause chronic illnesses and any age changes.

Shortage  $\text{CO}_2$  - a starting point for various disorders in an organism. If it proceeds during long time, it can be responsible for illnesses and ageing.

With the years at us physical activity that leads to reduction of carbonic gas in blood decreases and causes narrowing vessels. It's hard supply by oxygen of sites as a result of bad blood circulation in different sites of an organism. For example, if it occurs in a brain at the person arterial pressure if in heart the ischemia of heart and so on develops raises.

Normal maintenance  $\text{CO}_2$  in blood (6,5 %) provides a normal tone artery supports a good gleam of vessels.

It is possible to restore abilities of an organism to support optimum concentration  $\text{CO}_2$  in blood without additional, heavy physical activities Therefore is developed a method which is offered to you for maintenance of health of an organism.

Different methods of increase CO<sub>2</sub> are known: it is respiratory gymnastics jogi, run, a sauna, gymnastic of Butejko, Streljchenko and others.

Several years ago have invented the unique device for constant restoration CO<sub>2</sub>. Its name «Samozdraw» (to myself I do health).

The original idea was used and in the program of preparation of flights of cosmonauts.

This new device - for personal use of the house, also is a fast method of restoration of opportunities of an organism.

Everyone can use this device of the house, without supervision of the doctor.

To learn what from you percent of maintenance CO<sub>2</sub> in blood, the device includes the device, that to measure this parameter and if necessary to restore it.

Thus, the person, using the given device can restore the health, by means of not complex procedures.

How the device works?

Breathing by means of the device within 30-minutes it is possible to increase level CO<sub>2</sub> in an organism and it remains within several hours. It is a way to train the respiratory center to adapt to increased level CO<sub>2</sub>. After training several months, the respiratory center repeatedly adapts automatically to a normal level on 6.0 % CO<sub>2</sub>. Within these months, there is a reduction in depth of breath and ventilation of lungs, and level CO<sub>2</sub> in alveolar air in lungs gradually becomes normal (6.0 %). It in turn improves a metabolism in all cells of an organism, conducting to be reduced the help in all soft muscles (vessels, a bronchial tube ...) and increase of inviolability in general. The Human body could then to become healthy, and it is a lot of attributes of illness will disappear. That is more, cells become more proof to гипоксия because the stock of delivery of oxygen is created. The human body becomes more and more proof to an extreme pressure, and to physical, psychological and can prevent many illnesses, type of heart attack or an insult, a heart attack.

### Instructions on assembly capnicator

Capnicator serves as the simulator of physical activities of the person. Employment with it promote strengthening of bodies of breath, cardiovascular system, a gastroenteric path, the support-impellent device. Trainings with it aparate - preventive means for prevention of the diseases caused by a spasm of the main muscles and lack of oxygen in fabrics. The person, at breath through Capnicator forms an air mix of exhalation and external air.

### Rules for the training with capnicator.

1. During training You should breathe only Through a mouth.

2. Nose should be completely Blocked (fingers, Paper clips, wadded tampons, etc.).

3. Breathe easy, as usually, without efforts and delays, without efforts and delays. Duration of a session - 30 minutes. To breathe it is possible till a breakfast or for the night, in 2 hours after a supper

4. In the beginning it is not necessary to do procedure. With full duration (20 Minutes). Let the respiratory center will get used to loading. Begin with 7 minutes, and Step by step you will go To normal duration in 20 - 30 minutes.


5. Some people to be afraid of concept "To overtrain", on the some people stages. Not to be-the nature does not resolve search Infinite increase CO<sub>2</sub>, is a limit - 6,5% for Lives of alive fiber

6. Main parts of capnicator (see picture):

Base, 2 - Roof of trainer, 3 - Middle cell, 4 - Glass, 5 - Lid, 6 - Breathing hose, 7 - Mouthpiece .

### Sketch of assemblage

- The lid puts into the glass, that as opposite of mark on the glass will be the label 1 on the lid.
- Middle cell is putting in the base.

- The base is joining with roof. By this the projection on inside of rim of the roof (picture 2) to must fall in the groove on the base's rim (picture 5). The incidence of projection with groove is accompany with flick by turn of roof on the base.
- Middle cell to return such as her mark (picture 4) is inciting with the point of madding of trainer (picture 3) on the roof of trainer.
- The glass , which was gather with roof of the glass you must to put through the hole in the middle cell and to put his such , as marks on the glass is incites with the point of madding of trainer on the roof of trainer.
- The breathing hose is joining self her the end with projection of glass's lid, and second end - with bush of mouthpiece.

The training with capnicator consist of four stages, the month's duration each. On each stage capnicator must be reassembled.

#### 1-st stage

- Lid is inserted into glass so, that glass label matches with number 1 on the lid's top.
- Put breathing hose on lid. Put mouthpiece on breathing hose.
- It is all. Others parts (middle cell, roof, base) are unnecessary on 1-st stage!

#### 2-st stage (pic.5)


- The trainer is gathering according of sketch of assemblage.
- Middle cell to return such, as her mark is concur with the point II on the roof of the trainer 2.

#### 3-d stage (pic.6)

- The trainer is gathering according of sketch of assemblage.
- Middle cell to return such, as her mark is concur with the point III on the roof of the trainer 2.

#### 4-th stage

- Capnicator is assembled as on 3-d stage. It is poured 3 spoon of warm water into the glass. Lid of the glass is inserted such, as the mark of the glass will incite with on the label 1 on the lid of the glass (as on previous stages).
- Next 15 days lid is rotated on label 2.
- Next 15 days lid is rotated on label 3.
- Next 15 days lid is rotated on label 4.

## Measuring CO<sub>2</sub> in arterial blood

Yes, CO<sub>2</sub> is very important gas for the man's organism, but how it can be measured?

The units of measurement can be expressed in:

- by percentage of CO<sub>2</sub> in one unit of gas mixture. For example, atmosphere air contains 0,03% of CO<sub>2</sub>. Norm in alveolar air in lungs is 6%, in arterial blood - 6.5%, in cells 7%.
- by relative value, so 1 liter of arterial blood must contain 500 ml of CO<sub>2</sub> i.e. norm is 500 ml/L.
- millimoles per liter - relative quantity of molecules of CO<sub>2</sub> in 1 liter, example - 25mmol/L.
- partial pressure of CO<sub>2</sub> (PaCO<sub>2</sub>), example - 40 mm Hg.

The last two are traditionally getting as a result of laboratory testing. The blood sample is obtained by arterial puncture, usually in the wrist.

Normal blood gas results are as follows:

- Partial pressure of oxygen (PaO<sub>2</sub>): 75-100 mm Hg
- Partial pressure of carbon dioxide (PaCO<sub>2</sub>): 35-45 mm Hg
- Oxygen content (O<sub>2</sub>CT): 15-23%
- Oxygen saturation (SaO<sub>2</sub>): 94-100%
- Bicarbonate (HCO<sub>3</sub>): 22-26 mEq/liter
- PH: 7.35-7.45

Values that differ from those listed above may indicate respiratory, metabolic, or kidney disease.


*Despite the fact that PaCO<sub>2</sub> from 45 down to 35 mm Hg is considered to be normal, it is known that PaCO<sub>2</sub> reduction on 1 mm Hg -will change for the worse the blood supply on 3-4%. So, the reduction PaCO<sub>2</sub> by 10 mm Hg (to lower border of 35 mm Hg) will result to worsening of brain blood supply on 30-40%. Too bad!*

## Measurement of CO<sub>2</sub> with capnometer.

Blood gas analysis measures the partial pressures of oxygen and carbon dioxide in the blood, as well as oxygen content, oxygen saturation, bicarbonate content, and blood pH.

But we are interested in measuring CO<sub>2</sub> only. There is precise relation between the intensity of breathing and the level of CO<sub>2</sub> in arterial blood (the formulas are described in the physiology tutors). So, with the help of those formulas, the table has been made, using which you can determine your content of CO<sub>2</sub> by measuring your breath intensity.

Attention: *the measurement must be conducted in the state of full rest - best time is immediately after sleep or after not less than 30 min. of full immobility.*

### The construction of capnometer.

Its construction is very simple - it is a sack with valves (flappers). The foil volume of sack equals 13 litres. In one corner of sack there is a valve of lowering for air used for preliminary discharge of the air from the sack. Before the beginning of measurement, capnometer it is necessary to close the valve of lowering for air. After measuring procedure it is necessary to open the valve of lowering for air and to pit air from the sack.

For measurement the sack is inflated on the full with normal breathing. By the time of its filling (1-4 min) it is possible to calculate (with the help of table) intensity of breath and then concentration of CO<sub>2</sub>.

Breathing tube (hose) from capnicator is used, which must be connected to valves distributor in capnometer for the time of measurement. Valves distributor consist of two valves, one for inhaling external air and another for exhaling into the sack.

### The sequence of actions for measurement.

- insert the breathing tube (hose) into the aperture of the valves distributor.
- let an air from capnometer out via the hole in the corner of the sack.
- take the tube in the mouth, *squeeze the hole* in the corner of the sack with fingers, close the nose with fingers (or clips, cotton wool, etc.)
- note the time (with the second's precision) and start breathing into the sack .....

Attention: *breathe calmly as usual, but only through mouth.*

- then sack is filled, note the time. Calculate time interval (1 - 3.2 min.) and with the help of table define your value of breath per minute (VBM) and level of CO<sub>2</sub> in blood.

### After measurement:

- write data (time interval, percentage of CO<sub>2</sub>, date of measurement) in your private table.
- let an air from capnometer out, take the breathing tube out from the valves distributor.

*Remember to make measurements once or twice a month.*

Time Min. sec	VBM Litre / mia	CO <sub>2</sub> %
1m 00s	13.0	3.3
1m 10s	11.1	3.6
1m 20s	9.8	3.8
1m 30s	8.7	4.0
1m 40s	7.8	4.2
1m 50s	7.1	4.4
2m 00s	6.5	4.6
2m 10s	6.0	4.8
2m 20s	5.6	5.0
2m 30s	5.2	5.2
2m 40s	4.9	5.4
2m 50s	4.6	5.6
3m 00s	4.3	5.8
3m 10s	4.1	6.0
3m 20s	3.9	6.2

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## RHYSICAL TRAINER MIMIC (CAPNICATOR) Destis PASSPORT

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### 1. Appointment.

Trainer intends for realization of trainings the the organism of man , in the way creation in the organism of man the physiological condition, identi-Cal for condition , wich is arising in the resulte of physical exercises

Mimicking of physical exercises consist in increase of concentration carbon dioxide ( $CO_2$ ) in arterial blood. There is in process of recurrent respire through trainer , by this the concentration  $CO_2$  in the breath air will be increase. It depend on mood at trainer.

### 2. Complete set.

Base	1 piece
Roof of trainer	1 piece
Middle chamber	1 piece
Glass	1 piece
Lid	1 piece
Breathing hose	1 piece
Mouthpiece	1 piece
Clamp for nose	1 piece
Capnometer	2 piece
Passport and instruction on.	1 piece

### 3. Technics safe.

Over to keep the trainer in the package or in the paste box. Technical conditions on the trainer and trainer self have sanitary- epidemiological rules and normatives . All detales can be sterlisate\*of 3 % solution hydropririt on OCT -- 42-21-2-85 -- by immersion on 30 minutes with 0,5 % solution of domestics means by 18 -24° C temperature. The number of certificate 63. 01. 05. 961. П. 004040. 12.01 in 03.12.2001.

### 4. Evidence about receiving.

The trainer correspond to technical conditions. The date of output - December 2007.

### 5. Guarantee of maker.

The maker is guaranting the correction of trainer in time 18 monthes with making date or 15 monthes with sale's date.

### 6. The rule is producing of reclamation.

Reclamations can be prducing on conditions rule by under taking- making on adress : 443051 .Samara. F.niscyskaja street , 62-A , Scientific- making undertaking " Samozdrav". Tel. ( 846 ) 931-42-76.

### Attention !

Necessary all detales of trainer to wash with warm water with little dose of hold soda or detergent means before first use.

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